PACOTE CALIBRADOR ÍONPRO

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IONPRO FDS Calibrator Pack (ION PROLYTE Calibration Pack).

PRODUCT DESCRIPTION

The ION PROLYTE Calibration Pack is a set of calibrator solutions developed for use in ion-selective electrode (ISE) analyzers. The product consists of two standard solutions (Standard A and Standard B), used during the instrument's automatic calibration processes, ensuring correct electrode adjustment and the analytical performance of the system. The product is supplied ready to use and is intended exclusively for in vitro diagnostic use.

INTENDED USE

The ION PROLYTE Calibration Pack is intended for the automatic calibration of ion-selective electrodes for the quantitative determination of sodium (Na^+), potassium (K^+), chloride (Cl^-), ionized calcium (Ca^{2+}) and pH in compatible analyzers. Its use ensures the correct adjustment of the measuring system, contributing to the reliability of analytical results when used in accordance with these Instructions for Use.

This product is intended exclusively for in vitro diagnostic (IVD) use.

COMPATIBLE INSTRUMENTS

The product must be used exclusively on the following instruments:

Compatible instruments
XI-921
XI-931
XI-1021
IonPro

OPERATING PRINCIPLE

The system uses two reference calibrator solutions (Standard A and Standard B) containing known concentrations of the analytes. During automatic calibration, the ion-selective electrodes are sequentially exposed to the two standard solutions. The electrical potential difference obtained between the standards is used by the analyzer to calculate the slope and the reference point (offset) of each electrode, establishing the calibration curve required to determine the concentration of the analytes in patient samples.

This procedure ensures that the instrument operates within the analytical specifications established by the manufacturer.

COMPOSITION

Standard A

Component	Concentration
Sodium	< 180 mmol/L
Potassium	< 20 mmol/L
Chloride	< 100 mmol/L
Calcium	< 5 mmol/L
Preservative	< 1%

Standard B

Component	Concentration
Sodium	< 120 mmol/L
Potassium	< 20 mmol/L
Chloride	< 100 mmol/L
Calcium	< 5 mmol/L
Preservative	< 1%

PRESENTATION

Product	Presentation
ION PROLYTE Pack	Standard A + Standard B

STORAGE AND TRANSPORT CONDITIONS

Condition	Temperature range
Storage	15 °C to 30 °C
Transport	15 °C to 30 °C

Additional conditions:

- Keep protected from direct sunlight.
- Avoid exposure to excessive humidity.
- Do not freeze.
- Keep the package properly closed when not in use.

STABILITY AND SHELF LIFE

The product remains stable until the expiry date indicated on the package when stored under the specified conditions.

After opening and installation on the instrument, use within 60 days is recommended, always respecting the expiry date printed on the label.

Do not use expired reagents.

REAGENT PREPARATION

The reagent is supplied ready to use; no dilution or prior preparation is required.

PROCEDURE FOR USE

Pack installation

1. Check that the package is intact, correctly labelled and within its shelf life.
2. Carefully remove the grey silicone seal from the package, avoiding damage to the connector.
3. Place the Pack in the reagent compartment, according to the instrument manufacturer's instructions.
4. Connect the aspiration line, making sure it is fully seated, with no gaps that would allow air to enter.
5. Check all connections, making sure there are no leaks, obstructions, kinks or disconnections.
6. After installing the reagent, perform the Registration procedure on the instrument.
7. Perform a new automatic calibration.
8. Wait for the complete removal of air from the tubing.
9. Confirm that the calibration results meet the criteria established by the instrument before starting the analytical routine.

LIMITATIONS OF USE

- Product intended exclusively for the calibration of compatible instruments.
- Do not use on instruments not specified in these Instructions for Use.
- Do not reuse empty packages.
- Do not mix with other reagents.
- Do not transfer the contents to another container.

- Do not use reagents showing signs of leakage, contamination or visual alteration.

WARNINGS AND PRECAUTIONS

1. Product intended exclusively for in vitro diagnostic use.
2. Strictly follow the instructions described in this document.
3. Use personal protective equipment during handling.
4. Avoid contact with eyes, skin and mucous membranes.
5. In case of accidental contact, wash immediately with plenty of water.
6. The product may be harmful if swallowed or inhaled.
7. If swallowed, seek medical attention immediately.
8. Do not use damaged packages.
9. Instruments must be properly calibrated and subjected to periodic maintenance.
10. Disposal must comply with current environmental and health regulations.
11. Additional safety information is available in the Safety Data Sheet (SDS) at www.lscientific.com.br.

REFERENCE VALUES

Sample	Na^+ (mmol/L)	K^+ (mmol/L)	Ca^{2+} (mmol/L)	Cl^- (mmol/L)	pH
Serum / Plasma / Blood	136 - 145	3.5 - 5.1	1.12 - 1.32	97 - 111	7.35 - 7.45
Urine	40 - 220	25 - 125	Not applicable	110 - 250	Not applicable

DISPOSAL

Remaining reagent, used packages and other waste generated during use of the product must be disposed of in accordance with the institution's internal procedures and current health and environmental regulations.

Do not discharge the contents directly into sinks, sewage systems, soil, watercourses or ordinary waste without prior assessment of the applicable local regulations.

Packages containing product residues must be treated as chemical waste from healthcare services. Where the package, tubing or reagent has come into contact with biological samples, the material must be considered potentially contaminated and disposed of as healthcare waste, in accordance with the classification and procedures established by the institution.

Do not reuse, puncture, cut or use the empty package for other purposes.

QUALITY ASSURANCE

All lots of the ION PROLYTE Calibration Pack undergo internal Quality Control procedures before release for sale.

Product quality is guaranteed until the expiry date printed on the package, provided that the storage, transport and use conditions specified in these Instructions for Use are respected.

CONSUMER INFORMATION

LOT NUMBER, MANUFACTURING DATE AND EXPIRY DATE: SEE PRODUCT LABEL.

MANUFACTURER: LS CIENTIFICA LTDA - CNPJ: 03.611.091/0001-90
Rua Rio Parnaíba, 347 - Eldorado - CEP: 32371-380 - Contagem - MG -
Brazil
Phone: +55 (31) 3353-2173 - E-mail: lscientifica@lscientifica.com.br
Technical Manager: Wilson Adevaldo Inácio Silva - CRBM-MG: 15.974
Customer Service: +55 (31) 3353-2173

DISTRIBUTED BY: FDS TECNOLOGIA LTDA. CNPJ: 25.425.163/0001-47
Avenida Anápolis, 100, sala 1110, Bethaville - CEP: 06404-250 - Barueri/SP
- Brazil
Technical Manager: Jenifer Duarte Lima - CRF/SP: 73263 Anvisa: 8298271

Contact information (Customer Service):

Phone: +55-11-5200-0810
E-mail: contato@fdstec.com.br

PACK Family

Standard A - ANVISA notification No.: 80210830021
Standard B - ANVISA notification No.: 80210830023

Revision 01 - October 2026

MEANING OF THE SYMBOLS USED ON THE PRODUCT LABEL	
	Contains sufficient for <n> tests
	Use-by date
	Temperature limit (store at)
	Consult instructions for use
	Caution
	In vitro diagnostic medical device
	Biological risk
	Date of manufacture
	Batch code
	Catalogue number
	Do not re-use
	Manufacturer